



## Freestyle

- ✓ easy and fast assembly
- ✓ slim design: 9 cm overall height
- ✓ high stability

The comfortable and electrically adjustable Freestyle slatted frame consists of 42 wooden spring slats and has an overall height of 9 cm. The Freestyle can be used with all bed frames.



Back rest  
0 - 70°



Leg rest  
0 - 30°

### As an option

- wireless hand set
- mattress brackets
- steel side rails
- lifting pole
- simple extension to 210 or 220 cm length

### Standard features

- 42-slats comfort bed base
- max. patient weight 135 kg, safe working load 185 kg
- electrically operated back rest up to 70° and leg rest up to 30°
- overall height including motor: 9 cm
- 29 volt power supply for extra safety
- protection class IPX4



### Would you like to learn more?



Link to product page

Did we raise your interest?

For further information,  
please visit our website

[tecfor-care.com](https://www.tecfor-care.com)

We're happy to help you  
with your inquiry.

#### platform base

80 x 200 cm  
90 x 200 cm  
100 x 200 cm  
120 x 200 cm  
140 x 200 cm

#### outer dimensions

79 x 193,5 cm  
89 x 193,5 cm  
99 x 193,5 cm  
119 x 193,5 cm  
139 x 193,5 cm

1/2