



Care Bed Seracare

- ✓ central brake locking of the rear wheel
- ✓ preprogrammed comfort position
- ✓ individual designs of the bed ends
- ✓ auto regression of back rest
- ✓ many individual options

The “all-in-one care bed” combines all the requirements of a modern care bed in terms of comfort and ergonomics, and comes with a range of stylish designs. The Seracare is particularly suitable for use in care homes.



Back rest
0 – 70°



Leg rest
0 – 35°



Height
240 – 800 mm



Trendelenburg
included

Standard features

- 4-section steel slats platform base 90 x 200 cm with mattress retainer at head and leg rest
- electrically operated height adjustment 240 - 800 mm
- max. patient weight 190 kg, safe working load 225 kg
- electrically operated back rest up to 70° and leg rest up to 35°
- mechanical height adjustment of foot rest in 5 positions
- Trendelenburg / Anti-Trendelenburg function
- 4x castors 75 mm, 2x castors 100 mm, with central locking
- 35 volt power supply for extra safety
- 9V battery for emergency
- handset with single fault safety and locking function
- protection class IPX4
- auto regression of back rest up to 76 mm

As an option

- integrated bed extension up to 20 cm
- various decors
- various head- and footboard designs
- side rail options: split aluminium side rails, split steel side rails, continuous wooden side rails 95 mm



platform base
90 x 200 cm
100 x 200 cm


outer dimensions
105 x 209 cm
115 x 209 cm

Care Bed Seracare

Available decors

-  Altmuehl Beech
-  Lindberg Oak
-  More
-  Choco Bella Noce

Accessories for this bed



Would you like to learn more?



Link to product page

Did we raise your interest?

For further information,
please visit our website

tecfor-care.com

We're happy to help you
with your inquiry.