



Ecofit S

care bed
standard version



Ecofit S

The Ecofit S care bed has been specially designed for home care. Its high level of functionality and durability make it the ideal product in this sector.



Trendelenburg
optional



Back rest
0 - 70°



Hight
385 - 810 mm



Leg rest
0 - 30°


Standard features

- 4-section steel slats platform base with adjustable mattress retainers at head and leg rest
- electrically operated hight adjustment 385-810 mm
- max. patient weight 185 kg, safe working load 220 kg
- electrically operated back rest up to 70° and leg rest up to 30°
- mechanical height adjustment of foot rest in 3 positions
- 4x castors 100 mm, individually lockable
- wooden siderails (115 x 23 mm)
- 35 volt power supply for extra safety
- 9 volt battery for emergency
- handset with single fault safety and locking function
- protection class IPX 4
- patient lifting pole that can support up to 80 kg
- weight of each component is less than 20 kg

Accessories



different side rails



lots more accessories

- 20 cm bed extension
- handset with Trendelenburg/
Antitrendeleburg function
 - various decors
 - mattress levelling
- integrated bed extension

More upgrade and customisation options are available on the tecfor care upgrade-path.

As an option and resource numbers

model	platform base	outer dimensions
Ecofit S 80	80 x 200 cm	92 x 212,3 cm
Ecofit S	90 x 200 cm	102 x 210,9 cm
Ecofit S 100	100 x 200 cm	112 x 210,9 cm
Ecofit S 120	120 x 200 cm	132 x 212,8 cm

model	HMV-No.	Care-HMV-No.
Ecofit S	19.40.01.3268	50.45.01.1306
Ecofit S 100	19.40.01.3040	50.45.01.1079